# LUX-3D CAST

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# LUX 3D CAST

## DESCRIPTION

LUX 3D CAST is a light-curing, calcineable resin specially formulated for the impression of dental skeletons, crowns and bridges using the DLP AND MSLA additive manufacturing technique. This resin cures efficiently using a 385/405 nm UV-LED light source.

It burns out without leaving residues, making it the ideal choice for achieving precise and porosity-free castings. This feature ensures that the resulting dental pieces are of the highest quality and accuracy.

Once the printing process is completed, the pieces exhibit exceptional dimensional stability. This property is crucial for preserving the integrity of the most delicate and complex details, thus ensuring precise results and an optimal fit.

The LUX 3D CAST resin offers the possibility to continue the casting process using the CALIBRA 3D investment material, specially designed for this purpose.

## PROPERTIES

- Minimal shrinkage
- Burnout, leaving no residues.
- Precise reproduction of details.
- High surface hardness, providing dimensional stability.

## TECHNICAL DATA

• Viscosity:	<b>100 mPa.s</b>
• Flexural strength:	<b>60 MPa</b>
• Hardness:	<b>75 Shore D</b>
• Layer thickness:	<b>From 30µm to 100µm.</b>
• Colour:	<b>Red</b>
• Shrinkage:	<b>&lt;0.5%</b>
• Ash content:	<b>&lt;0.1%</b>

## COMPOSITION

- Monomers
- Photo initiators
- Additives

## INSTRUCTIONS FOR USE

### WORKING PARAMETERS

#### PRINT PREPARATION

- Design the dental restoration according to your CAD software instructions, considering the minimum thickness requirements.
- Once the design is completed, import the STL file into your slicing program. Place the file and corresponding support structures as per your printer manufacturer's instructions.
- Enter the working parameters based on your printer model. You can find tested and compatible 3D printing systems on our website.

#### [PRINTER CONFIGURATION](#)

*If you have another printer, please contact us. We can help you with the parameterisation of your printer.*


#### MODEL PRINTING

- Make sure that the temperature of the resin and the room is within the range of 18 to 28 °C.
- Before using LUX 3D cast, shake the container for approximately 2 minutes to ensure a homogeneous mixture of the components. (Fig.1)




- To prevent possible deformations and printing errors, make sure all working instruments are clean.
- Pour the resin into the vat, and before starting the printing, wait a few minutes for the bubbles to rise to the surface. (Fig.2)

- After the printing process, we recommend letting the objects drain for 10 minutes.
- Return the excess material to a separate container, do not mix it with unused resin. Make sure the material does not contain impurities or polymerized residues using a paper filter.

## CLEANING

• The properties of the finished product depend, among other factors, on post-processing.

• Cleaning should be done in two steps:

### Pre-cleaning:

Subject the printed objects to a pre-cleaning process for 2-3 minutes in a reusable isopropyl alcohol bath. **(REF. LUX 3D CLEAN)**. (Fig.3)


Note: The effectiveness of the cleaning bath diminishes with each use. Therefore, replace the bath if the cleaning effectiveness is reduced. Thoroughly clean the cavities and hollows, if necessary, using compressed air.

### Main cleaning:

Next, clean the printed objects for **6 minutes** in a clean isopropyl alcohol ultrasonic bath (purity >98%). **(REF. LUX 3D CLEAN)**. (Fig.4)

Finally, thoroughly dry them using compressed air.

## RE-HARDENING - POST CURING


• Post-curing is performed using post-curing machines.

• The post-curing times validated with our machines are as follows:

Machine	Programme	Time
OTOFLASH	2x2000 Flashes	
BB CURE	Clear	20 minutes

• Any other machine to be used requires its own tests to validate the process.

## FINAL FINISHING

• Carefully remove and cut the support structures using a non-water-cooled rotary tool or cutting pliers without applying force. (Fig.5)

• Use a suction device. Thoroughly remove any remaining resin dust with compressed air.

## OBSERVATIONS, PREVENTIVE MEASURES

Please refer to the safety data sheet for LUX 3D CAST resin! (Available on the web.)

LUX 3D CAST contains acrylates. LUX 3D CAST should not be used if you have hypersensitivities (allergies) to these ingredients.

Use appropriate personal protective equipment (gloves and protective goggles) during use. Avoid direct contact with the liquid material and its components before curing.

## STORAGE


Store at a temperature between 15°C and 28°C.

The batch number and expiration date are marked on each container.

Close the jar immediately after each application.

The material cures when exposed to light.

Do not use the product after the expiration date has passed.

## WASTE MANAGEMENT

The polymerized material can be safely disposed of along with household waste. However, it is important to note that unused resin and isopropyl alcohol used for cleaning, containing resin residues, should not be disposed of with regular household garbage. These waste materials should not be poured down the drain but should be disposed of according to applicable local regulations for safe and proper disposal. It is important to follow regulations and recommended practices to protect the environment and public health.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

PARA USO DENTAL EXCLUSIVAMENTE.

FOR DENTAL USE ONLY.

Polígono Empordà Internacional  
17469 VILAMALLA, GIRONA • SPAIN  
Tel.: ( 34 ) 972 52 61 69  
protechno@protechno.com  
**www.protechno.com**