IMPRESIÓN 3D • 3D PRINTING

LUX-3D MODEL 4K

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LUX 3D MODEL 4K

DESCRIPTION

LUX 3D MODEL 4K is a resin specifically designed for 3D printing of orthodontic dental models and prostheses using additive manufacturing techniques such as DLP and MSLA. This resin cures with a 385/405 nm UV-LED light source, ensuring precise and detailed printing. Furthermore, LUX 3D MODEL 4K is resistant to shrinkage and offers high surface hardness and scratch resistance, providing dimensional stability and durability to dental models.

LUX 3D MODEL 4K features an opaque color that allows for precise optical recognition of preparation boundaries, facilitating the design and production process of high-quality dental models.

It is important to note that this resin is safe for use as it is completely free from the hazards of TPO (Diphenyl (2,4,6-trimethylbenzoyl) phosphine oxide (CAS 75980-60-8)).

LUX 3D MODEL 4K offers a range of unique advantages for model printing, including its ability to withstand thermoforming processes. This means the resin can withstand high temperatures and pressures during the shaping process, resulting in more accurate and higher-quality models.

PROPERTIES

- Minimal shrinkage < 0.5%
- Precise reproduction of details.
- High surface hardness and scratch resistance, providing dimensional stability.
- Opaque color adjustment for optical recognition of surface structure, preparation boundaries, etc.
- Suitable for thermoforming processes.

TECHNICAL DATA

- Viscosity: **1200 mPa.s**
- Flexural strength: **80 MPa**
- Hardness: **80 Shore D**
- Elastic modulus: **1500 MPa**
- Colors: apricot, white, and anthracite.

The values were obtained using Lux 3D printer and standard printing parameters.

COMPOSITION

- Monomers
- Photo initiators
- Additives

INSTRUCTIONS FOR USE

WORKING PARAMETERS

PRINT PREPARATION


- Design the model according to the instructions provided by your CAD software, considering the minimum wall thickness requirements. A minimum wall thickness of 3 mm is recommended for hollow models.
- Once the design is complete, import the STL file into your slicing program. Position the file and add the necessary support structures as per the instructions provided by your printer manufacturer.
- Enter the working parameters specific to your printer model. You can find the tested and compatible 3D printing systems on our website.

[PRINTER CONFIGURATION](#)

If you have another printer, please contact us. We can help you with the parameterisation of your printer.

MODEL PRINTING


- Make sure that the temperature of the resin and the room is within the range of 22 to 28 °C.
- Before using LUX 3D MODEL 4K, shake the container for approximately 2 minutes to ensure a homogeneous mixture of the components. (Fig.1)




- To prevent possible deformations and printing errors, make sure all working instruments are clean.

- Pour the resin into the vat, and before starting the printing, wait a few minutes for the bubbles to rise to the surface. (Fig.2)

- After the printing process, we recommend letting the objects drain for 10 minutes.
- Return the excess material to a separate container, do not mix it with unused resin. Make sure the material does not contain impurities or polymerized residues using a paper filter.

CLEANING

• The properties of the finished product depend, among other factors, on post-processing.

• Cleaning should be done in two steps:

Pre-cleaning:

Subject the printed objects to a pre-cleaning process for 2-3 minutes in a reusable isopropyl alcohol bath. **(REF. LUX 3D CLEAN)**. (Fig.3)

Note: The effectiveness of the cleaning bath diminishes with each use. Therefore, replace the bath if the cleaning effectiveness is reduced. Thoroughly clean the cavities and hollows, if necessary, using compressed air.

Main cleaning:

Next, clean the printed objects for **6 minutes** in a clean isopropyl alcohol ultrasonic bath (purity >98%). **(REF. LUX 3D CLEAN)**. (Fig.4)

Finally, thoroughly dry them using compressed air.

RE-HARDENING - POST CURING



- Post-curing is performed using post-curing machines.
- The post-curing times validated with our machines are as follows:

Machine	Programme	Time
OTOFLASH	6000 Flashes (with Nitrogen)	
BB CURE	Clear	20 minutes

- Any other machine to be used requires its own tests to validate the process.

FINAL FINISHING

• Carefully remove and cut the support structures using a non-water-cooled rotary tool or cutting pliers without applying force. (Fig.5)

• Use a suction device. Thoroughly remove any remaining resin dust with compressed air.

OBSERVATIONS, PREVENTIVE MEASURES

Please refer to the safety data sheet for LUX 3D MODEL 4K resin! (Available on the web.)

LUX 3D MODEL 4K contains acrylates. LUX 3D MODEL 4K should not be used if you have hypersensitivities (allergies) to these ingredients.

Use appropriate personal protective equipment (gloves and protective goggles) during use. Avoid direct contact with the liquid material and its components before curing.

STORAGE


Store at a temperature between 15°C and 28°C.

The batch number and expiration date are marked on each container.

Close the jar immediately after each application.

The material cures when exposed to light.

Do not use the product after the expiration date has passed.

WASTE MANAGEMENT

The polymerized material can be safely disposed of along with household waste. However, it is important to note that unused resin and isopropyl alcohol used for cleaning, containing resin residues, should not be disposed of with regular household garbage. These waste materials should not be poured down the drain but should be disposed of according to applicable local regulations for safe and proper disposal. It is important to follow regulations and recommended practices to protect the environment and public health.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

PARA USO DENTAL EXCLUSIVAMENTE.

FOR DENTAL USE ONLY.

Polígono Empordà Internacional
17469 VILAMALLA, GIRONA • SPAIN
Tel.: (34) 972 52 61 69
protechno@protechno.com
www.protechno.com